

# Side\_Pop\_Protocol

Callum Malcolm

2023-05-18

## Protocol

**Note:** Seed cells 1-2 days before so they are in log phase growth

**Note:** Preheat thermomixer to 37°C and take out Verapamil and Hoechst from -20°C before starting

1. Suspend cells and count
2. Dilute cells to  $1 \times 10^6$  cells/mL  
**Note:** Determine correct dilution media (FACS buffer, Media, ect.)  
**Ex:** For N2 use media as FACS buffer will reduce viability
3. Distribute 1ml of cell suspension to eppendorfs
4. Add X  $\mu$ L verapamil to control samples and X  $\mu$ L PBS to test samples  
**Note:** Verapamil volume depends on cell line; cell lines with larger morphology could require more volume. Important to optimise this aspect of the assay before proceeding.  
**Ex:** 10  $\mu$ L Verapamil is optimal for N2 PDX
5. Incubate samples for 20 min at 37°C shaking at 600rpm
6. Add Hoechst to samples  
**Note:** Be **VERY** precise with pipetting (ie. remove droplets from tip). This is the most important optimization step. Test with a range of concentrations with corresponding Verapamil controls to determine optimal volumes for both reagents.  
**Note:** Protect Hoechst from light (ie. cover shaker with tinfoil, try to operate in the dark)  
**Ex:** 2.5  $\mu$ L Hoescht is optimal for N2 PDX
7. Incubate samples in thermomixer for 2hrs at 37°C and 600rpm
8. Put samples on ice and spin down at 250 x g at 4°C  
**Note:** Set centrifuge to 4°C 30min before the end of Hoechst incubation, collect ice and label FACS tubes
9. Resuspend in ice cold dilution media  
**Note:** Use ~300  $\mu$ L
10. Read samples on Fortessa II

## Reagents

Reagent	Product Code	[Stock]	[Assay]	Volume of Stock Used
Hoescht	ab228551	1mg/ml	2-5 $\mu\text{g}/\text{mL}$	2-5 $\mu\text{L}/\text{mL}$
Verapamil	V4629	10mM	100 $\mu\text{M}$	10 $\mu\text{L}/\text{mL}$

## Cytometer Optics

- Email flow facilites about changing the filters before starting
- UV Laser (Fortessa II) excites at 350nm
- Dichroic mirror LP600 splits beam:
  - Hoechst Blue: 450/50 BP
  - Hoechst Red: 610/20 BP